

Order form

Date:

Client's data:

Rf.user*:

*You are required **NOT** to indicate personal and sensitive data of the user in accordance with General Data Protection Regulation EU 2016/679.

CONFIGURATION

Standing frame in prone and supine position

01 Models

PRONE SUPINE



P1
P2
P3

S1
S2
S3



02 Basic configuration

4 frame extensions with braked swivel castors ✓
high adjustment piston ✓
tilting piston ✓
2 chest support units ✓
2 positioning belts ✓
couple of knee cups ✓
couple of individual footrests M1/M2/M3 ✓
tray ✓

PRONE

SUPINE

4 frame extensions with braked swivel castors ✓
high adjustment piston ✓
tilting piston ✓
2 chest support units ✓
2 positioning belts ✓
couple of knee cups ✓
couple of individual footrests M1/M2/M3 ✓
tray ✓
headrest ✓

02 Accessories

WALLY-KIT-01 Couple of footrests M1=180
 WALLY-KIT-02 Couple of footrests M2=215
 WALLY-KIT-03 Couple of footrests M3=240
 WALLY-KIT-04 Headrest

WALLY-KIT-05 Headrest for prone position M3
 WALLY-KIT-06 Headrest for prone position M2
 WALLY-KIT-07 Headrest for prone position M1
 WALLY-KIT-08 Support tray

04 NOTES

For better device configuration, please indicate weight, height and age of the user.
Describe the user indicating diagnosis, posture characteristics or other relevant details.

① User's weight **kg** _____
② User's height **cm** _____
③ User's age **yrs** _____

(filled by healthcare professionals)

signature _____

TECHNICAL DATA		M1	M2	M3
User's age	yrs	1-5	5-13	10-14
Maximum user's height	cm	50-95	90-120	100-140
Maximum user's weight	kg	40	40	50
Device weight	kg	23	26	32
High adjustment piston	N	900	1100	1100
Tilting piston	N	600	600	600
Total length	cm	77-84	77-90	78-94
Total width	cm	66-68	68-74	76-79
Tilting block height (from ground)	cm	40 - 90	40 - 90	40 - 90
Tilting rage of orthosis	cm	90°	90°	90°

MOD.ORD.WALLY.eng.rev.01